


QUICK
master



Facebow

**INSTRUCTION
MANUAL**



FACEBOW INSTRUCTIONS



Placement and stabilization of the bitefork.





Placement and locking of facebow.




Placement of nasal guide and blocking.



Locking of rigid connection between facebow and bitefork.



Facebow in place.



With the Quick Facebow, "Put your patient in the articulator."

WWW.QUICKARTICULATOR.COM

SINTEC, INC.

350 Concord Lane • P.O. Box 129 • E. Wakefield, NH 03830
800-233-2162 • 603-522-8624